



**TIME MEDICAL**  
SYSTEMS



DISTRIBUTED BY  
**EURONOX**  
MEDICAL GROUP




**neona**  
neonatal MRI system

# neona

neonatal MRI system

The NEONA has most updated cryo-free magnet which does not require quench pipe. It also has compact footprint, light-weighted magnet and is specifically designed for installation within Neonatal Intensive Care Unit (NICU). It has incorporated with fast scanning protocols with reduced acoustic noise. The strong gradient ensures the high image quality in terms of high spatial resolution and signal-to-noise ratio under short acquisition time.

US FDA 510(k) is cleared.



# Features and Benefits



**1.5T Cryo-free**  
Superconducting Magnet



**No Quench Pipe**  
required



**Safety**  
Maintain a safe Specific  
Absorption Rate (SAR)



**Improve**  
image quality for diagnostic  
purpose on neonates



**Light weight & Compact**  
for installation on upper floors



Design of  
**Ultra-small FOV**  
RF receiver coils



**Specialized pulse sequences**  
for neonatal scans



**Noise reduction**  
Optimize gradient acoustic  
noise for best patient comfort

DISTRIBUTED BY



FOR ANY ENQUIRIES CONTACT US

---

020 8146 6839 | +44 75 2990 2750

[info@euronoxmedical.com](mailto:info@euronoxmedical.com)

[www.euronoxmedical.com](http://www.euronoxmedical.com)

35 Berkeley Square, Mayfair, London, W1J 5BF